



### Foot for Total Knee Replacement

- This new design offers the same secure support as the stainless steel model without the concerns of radiopacity
- The surgeon can replace the knee while the patient's foot is held firmly in place
- · Disposable conductive, non-latex foot pad included
- · Now available in radiolucent, clear Lexan®
- · Patient's well leg lies flat

#800-0011 Foot, Lexan #800-0062 Foot, Stainless Steel #508-0111 Disposable Foot Pads, 12/cs

#800-0134 Socket XL





#### Arthroscopic Leg Holder

- VELCRO® brand fastener strap wraps completely over leg holder for total stability
- · Pressure relief foam pad included
- · 5/8" (2cm) mounting post

Extremely versatile, easy-to-use and completely rigid. Spring-loaded braces slide into place and snap lock.

#800-0023 Arthroscopic Leg Holder #508-0103 Disposable Premium Pads, 12/cs





#### Grip Leg Holder

- Dual independent braces freely rotate 360° and spread horizontally from 3"-11"W (8cm-28cm) to fit any size leg at any angle
- Twist one brace side down, put the patient's leg on the leg holder, return the brace side back up, slide brace horizontally for a snug fit, and lock in place with convenient twist handle

Use also for ACL procedures. Stainless steel. 12"L x  $1\frac{1}{2}$ "D (31cm x 4cm) leg support x 8"H (20cm) braces, 16"H x  $\frac{5}{8}$ "D (41cm x 1.6cm) mounting post.

#800-0065 Grip Leg Holder



# TKR for Total Knee Replacement





## TKR Positioner for Total Knee Replacement

- · Track Locking Bracket securely locks the Base Plate to the table with two thumb screws
- · Locks in any position with a simple turn of the wrist
- · Smooth sliding motion for effortless leg positioning
- · All components are autoclavable
- · CamLoc Clamp sold separately (see below)

Ball joint allows the surgeon to position the operative leg in extension, flexion, tilt and rotation.

#800-0141 TKR Positioner, 24" (61cm)

#800-0171 TKR Positioner, 30" (76cm)

#800-0146 Disposable TKR Boot Liners, Sterile, 12/cs

Required Accessory: Choose One

CamLoc Clamp

#800-0271 US #800-0273 UK

#800-0272 EU

